

Whole genome assembly using long read sequences

Heidi Tschanz-Lischer, Interfaculty Bioinformatics Unit (IBU)

11.06.2019

Why do we need whole genome sequences?

- Sequencing the genome is an important step towards **understanding it**
- Better **understand variations** within and between species
- Makes it easier to study
 - Cause of diseases
 - Morphological variation
 - Environmental adaptation
 - Genomic basis for evolutionary speciation
 - Gene expression divergence
 - Epigenetic modifications
 - ...
- **Reduces the costs** of future sequencing projects
 - Lower coverage required → population genomic or genome wide association studies
 - Simplify bioinformatics analyses





Genome sequencing costs dropped

Next-Generation Sequencing

- Highly parallelized sequencing
- Reduces cost and time

→ **main challenge** in genome sequencing has shifted from data generation to **reconstruction of genomes**







Sequencing by synthesis (termination)

- 50-300 bp length (PE)
- Error rate: 0.1%
- Less susceptible to homopolymer errors
- Under-representation in AT-rich and GC-rich regions
- Tendency towards substitution errors
- High throughput: 15-3000 Gb


Single-molecule real-time (SMRT), PacBio



- Slowed down DNA polymerase
- Measure colored light emission

Single-molecule real-time (SMRT), PacBio



- Average read length (Sequel II): ~50 kb (up to reads of 175 kb length)
- Error rate: 10-15%
- Bias towards indel errors
- Errors are supposed to be random
→ high coverage can overcome high error rate
- Higher cost
- Unique circular template
 - allows each template to be sequenced multiple times
 - create a circular consensus sequence (HiFi reads)
 - reduces error rate to <1%
(10 passes: ~99.8% accuracy)



Oxford Nanopore sequencing

Nano-scale pores

- Drag DNA/RNA through the pore
- Electrical current flows through the hole
→ disruption is measured if DNA pass through
- Shifts in voltage
→ characteristic for particular sequence in pore (k-mer)
→ Modifications to the primary DNA or RNA (e.g. cytosine methylation) can be detected directly
- Problem:
Multiple bases influence the current passing through the pore




Oxford Nanopore sequencing



MinION: a portable sequencing device


- 512 channels
- Small (USB based device), runs off a personal computer
- real-time sequencing data (~150bp per seconds; 500bps fast mode)
→ no fixed run time: *Run until...* sufficient data (15-30 Gb)
- Long reads >6 kb (up to 2mb) → DNA molecule length dependent
- Error rate: 5-12% (dominated by indels, homopolymers)

R9: May 2016 Onwards



PromethION: small benchtop system

- contains docking for 48 flow cells, each with 3000 nanopores
→ Total: 144,000 nanopores (< 4.8-8.6 Tb)



VoITRAX v2: programmable, portable device for automatic sample and library preparation




De novo genome assembly

De novo assembly:

- Reconstructing the original DNA sequence from fragmented reads
 - is like a big and complicated **jigsaw puzzle**
 - Millions of small pieces
 - Missing pieces
 - Some pieces have mistakes (sequencing errors)
 - Polymorphisms (diploid)
 - Long repetitive parts
- Different algorithms and assemblers developed (e.g. ABySS, SOAPdenovo2, ALLPATHS-LG, Falcon, Canu, ...)



De novo assembly



What is needed for a good assembly?

- Low heterozygosity DNA
- High coverage
- High read lengths
- Good read quality

Current sequencing technologies do not have all


- **Illumina:** good quality reads, but short
- **PacBio / Nanopore:** very long reads, but lower quality

→ Genome assembly is still a difficult problem and requires high computational resources

Choosing assembly strategy

The choice of algorithms depends on

- how much long reads (PacBio/Nanopore) can be obtained
- how much short read data are available




Long read de novo assemblers

HGAP: Hierarchical Genome Assembly Process

(<https://github.com/PacificBiosciences/Bioinformatics-Training/wiki/HGAP>)

- PacBio, included in the SMART Analysis software (GUI based)
- developed to allow the complete and accurate assembly of bacterial sized genomes (<100 Mb)
- 3 step process
 - **Preadassembly:**
generate long and highly accurate reads
 - **Assembly:**
Overlap-layout-consensus (OLC)
 - **Consensus polishing:**
reduce remaining Indels and SNP errors (Quiver)



Long read de novo assemblers

Shasta (<https://github.com/chanzuckerberg/shasta>)

- Nanopore reads
- Very fast and simple to use
- Default parameters optimized for coverage ~60x
- Output assembly in FASTA and GFA 1.0 (assembly graphs)
- RAM requirements: around 5-8 bytes per input base
- Early indications are that Shasta accuracy is at least comparable to alternative assemblers
- Designed for *de novo* assembly of human genomes
 - RAM: 1 TB (60x coverage)
 - Runtime: 6 hours (128 virtual CPUs)

Long read de novo assemblers

Miniasm (<https://github.com/lh3/miniasm>)

- Very fast OLC-based (overlap-layout-consensus) de novo assembler for noisy long reads
- Outputs only assembly graphs (GFA format) → no consensus calling

Unicycler (<https://github.com/rrwick/Unicycler>)

- Illumina-only (SPAdes optimizer), hybrid or long-read-only (miniasm + Racon) assemblies
- Can cope with very repetitive genomes
- Not especially fast, but circularizes genomes without a separate tool (e.g. Circlator)

Both: → for bacterial genomes
→ PacBio or Nanopore reads
→ easy and straight forward to use



Long read de novo assemblers

Flye (<https://github.com/fenderglass/Flye>)

- PacBio or Nanopore
- from small bacterial projects to large mammalian-scale assemblies
 - E. coli (4.6 Mb) 50x PacBio: 2h CPU time, 2 Gb RAM
 - Human (2.9 Gb) 30x PacBio: 900h CPU time, 300 Gb RAM
 - Human (2.9 Gb) 35x Nanopore: 5000h CPU time, 600 Gb RAM
- complete pipeline: raw reads → polished contigs
- Include special mode for metagenome assembly
- Easy and straight forward to use

Long read de novo assemblers




FALCON / pb-assembly (<https://github.com/PacificBiosciences/pb-assembly>)

- PacBio Assembly tool Suite
- diploid-aware assembler → follows HGAP
- optimized for large genome assembly
- >30-50x per haplotype (highly heterozygous diploid → require the double)
- extensive configuration file required
 - not easy to understand parameters
 - A few example files → can be used as a basis for modification

FALCON-Unzip:

- phase the genome and perform phased-polishing with Arrow
- partially-phased primary contigs and fully-phased haplotigs (haplotypes)



Long read de novo assemblers

Canu (fork of Celera Assembler; <https://canu.readthedocs.io/en/latest/index.html>)

- PacBio or Nanopore
- 3 phases: correction → trimming (get high-quality sequences) → assembly
- follows the hierarchical genome assembly process (HGAP)
- >30-60x
- automatically takes full advantage of grid systems (cluster) → submitting itself for execution
- consensus sequences:
 - >99% identity for PacBio
 - >98% identity for Nanopore (accuracy varies depending on pore and basecaller version)
- Easy and straight forward to use
 - Good manual with recommendations for parameter values (PacBio, Nanopore, low coverage data)

Post-assembly correction

→ improves quality and removes errors

Polish assembly with **long reads**

- Nanopore → **Nanopolish**:
 - calculates an improved consensus sequence
 - nanopolish call-methylation: predict methylated genomic bases
 - nanopolish variants: detect SNPs and indels
- PacBio → **Arrow** (former Quiver):
 - Get improved consensus → based on a hidden Markov model approach
 - get variant calls

Polish assembly with **Illumina reads**:

- **Pilon**:
 - Automatically improve draft assemblies (SNPs, small/large indels, gap filling, local misassemblies)
 - Find variations, including large event detection



Bacteria

Prokka

- rapid prokaryotic genome annotation
- quickly annotate bacterial, archaeal and viral genomes
- Outputs standard-compliant files

RAST

- Rapid Annotation using Subsystem Technology
- fully-automated service for annotating complete or nearly complete bacterial and archaeal genomes
- Webservice (<http://rast.theseed.org/FIG/rast.cgi>)

Eukaryotes

- **RepeatMasker**: screens genome for interspersed repeats and low complexity DNA sequences
- **MAKER**:
 - Genome annotation pipeline → allow smaller projects to independently annotate their genomes
 - identifies repeats
 - aligns ESTs and proteins to a genome
 - produces ab-initio gene predictions
 - especially useful for projects with minimal bioinformatics expertise and computer resources
- **PASA** (Program to Assemble Spliced Alignments):
 - exploits spliced alignments of transcripts to automatically model gene structures and splice variations
- **Augustus**:
 - find genes and their structures
 - can be used as an ab initio program → bases its prediction purely on the sequence.
 - also incorporate hints from extrinsic sources (e.g.: EST, MS/MS, protein alignments, ...)


Example



Cichlid genome assembly

- Cichlid lineage of the East African Great Lakes is famous for forming large adaptive radiations in exceptionally short time
- Recently, several cichlid genomes were published (*Brawand et al. 2014, Feulner et al. 2018*)

→ Improved genome will be useful for future studies of adaptation and speciation of Lake Victoria cichlids



Cichlid genome resources

Pundamilia nyererei



PunNye1.0 (Broad Institute):

- Nb scaffolds: 7,236
- N50: 2.5 Mb
- Total length: 830.1 Mb
- Total length: 698.8 Mb (without N)

- 126x Illumina read
- ALLPATHS-LG



PunNye2.0 (Feulner *et al.*):

- Nb scaffolds: 6,876
- N50: 29.8 Mb
- Total length: 856.2 Mb
- Total length: 698.8 Mb (without N)

- Linkage map - 1,597 SNP markers
- ALLMAPS

Raw data


PacBio (Sequel)

- Nb reads: 4,020,155
- Min length: 50 bp
- Max length: 143,514 bp
- Mean length: 10,538 bp
- Total length: 42.35 Gb → estimated coverage 42.7x



Illumina reads

- 4 closely related samples (380bp insertion):
 - Nb reads: 520,955,224
 - Total length: 78.14 Gb → estimated coverage ~78.7x (each sample 15-20x)
- SRA samples (used in original PunNye1.0 assembly):
 - 3 kb libraries: 709,783,284 (72.2x coverage)
 - 6-14 kb libraries: 721,087,418 (51.2x coverage)
 - 40 kb FOSILLs4: 36,341,216 (3.7x coverage)




Pipeline – genome assembly

Nb reads: 4,020,155
Mean length: 10,538 bp
Total length: 42.36 Gb (~42.7x)

Nb reads: 2,228,798
Mean length: 13,505 bp
Total length: 30.10 Gb (~30.3x)

Nb contigs: 6,732
N50: 0.920 Mb
Total length: 1.111 Gb

Nb contigs: 6,584
N50: 0.908 Mb
Total length: 1.083 Gb




Pipeline – genome assembly

Nb contigs: 6,584
N50: 0.908 Mb
Total length: 1.083 Gb

Nb scaffolds: 6,250
N50: 1.212 Mb
Total length: 1.085 Gb

Nb scaffolds: 5,753
N50: 27.607 Mb
Total length: 1.130 Gb
(without N: 1.083 Gb)






Nb contigs: 3,694
N50: 1.023 Mb
Total length: 1.033 Gb


Nb scaffolds: 3,335
N50: 1.324 Mb
Total length: 1.036 Gb

Nb scaffolds: 2,832
N50: 29.594 Mb
Total length: 1.074 Gb
(without N: 1.033 Gb)

Assembly comparison




Assembly comparison



Assembly comparison

	Stats	Pnyererei2	canu	falcon
QUAST				
nb scaffolds		6'876	5'753	2'832
N50		29'830'996	27'606'787	29'593'783
NG50		27'967'145	28'136'703	31'084'370
max length		60'199'168	62'610'257	55'908'494
total length		856'242'559	1'130'373'166	1'073'822'959
total length without N		698'778'000	1'083'432'443	1'032'560'317
N's per 100kb		18'390	4'153	3'843
BUSCO				
complete		2'498	2'485	2'527
complete single copy		2'469	2'255	2'391
complete duplicated		29	230	136
fragmented		58	38	37
missing		30	63	22

BUSCO Assessment Results



Conclusion

- Long read sequences are **important** for **high quality draft genomes**
- **>50x coverage** is required for long reads only assemblies
 - still expensive for large genomes, but prices will come down even more in near future
- **No single best assembly strategy/program**, depends on
 - Input data (quality, coverage)
 - Species (heterozygosity, complexity)
- Assembly evaluations are not straight forward
 - Longer assemblies (higher N50/NG50) are not always the best assemblies
 - Always use a combination of metrics
 - Only a few tools work without a known reference
(e.g.: BUSCO (Waterhouse et al. 2017), QUAST (Gurevich et al. 2013),
ALE (Clark et al. 2013), REAPR (Hunt et al. 2013))